Fosfatasa Acida

Total y Prostática cinética

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Fosfatasa Acida

Total y Prostática cinética

Colorimetric method for the determination of acid phosphatase and prostatic acid phosphatase

SUMMARY

Acid phosphatases are present in almost all body tissues, and the highest levels are found specially high in prostate, stomach, liver, muscle, spleen, erythrocytes, and platelets. The various isoenzymes differentiate each other by their optimal pH, molecular weight, activators and inhibitors.

Prostatic acid phosphatase (PAP) is a valuable method for the diagnosis of prostate cancer in an early stage, one of the neoplastic forms of higher levels of morbidity. PAP enzymatic activity using α -naphthyl phosphate as substrate, has shown adequate sensitivity and specificity for the diagnosis of prostate cancer.

Increases in the total acid phosphatase (ACP) activity have been detected in some hematological diseases (myelocytic leukemia, idiopathic thrombocytopenia) and bone diseases (Paget disease, bone carcinoma), as well as some types of cancer, hepatic diseases (hepatitis, obstructive jaundice), etc.

PRINCIPLE

The acid phosphatase (ACP EC 3.1.3.2.) hydrolyzes the α -naphthyl phosphate at pH 5.2 releasing phosphate and α -naphthol. Naphtol reacts with a diazo reagent (Fast Red TR) yielding a yellow pigment. Thus, the absorbance measured at 405 nm increases and it is proportional to the total acid phosphatase (ACP) levels in the sample. When the activity is measured in presence of tartrate, the prostatic isoenzyme activity is inhibited. The difference among the diverse activities of total acid phosphatase (ACP) and the one tartrate resistant (Non-Prostatic Acid Phosphatase (ACP-NP) corresponds to prostatic fraction (PAP).

PROVIDED REAGENTS

- A. Reagent A:** 0.1 mol/l citrate buffer, pH 5.2.
B. Reagent B: vials containing 10 mmol/l α -naphthyl phosphate and 1.5 mmol/l Fast Red TR.
C. Reagent C: 135 mmol/l tartrate solution in 100 mmol/l citrate buffer, pH 5.2.
D. Reagent D: 0.7 mol/l acetic acid solution.

NON-PROVIDED REAGENTS

Saline solution (9 g/l sodium chloride).

INSTRUCTIONS FOR USE

Total Acid Phosphatase (ACP)

Reagent B; preparation: dissolve with Reagent A as indicated volume in the label. Label ACP.

Non-Prostatic Acid Phosphatase (ACP-NP)

Reagent B; preparation: dissolve with Reagent C volume as indicated in the label. Label ACP-NP.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until their expiration date stated on the box. Room temperature exposure for extended periods of time may deteriorate Reagent B.

Reconstituted Reagent B: stable for up to 3 days at 2-10°C or for up to 24 hours at room temperature (< 25°C).

INSTABILITY OR DETERIORATION OF REAGENTS

If O.D. readings of Reagent B at 405 nm are above 0.250, it is an indicator of Reagent deterioration.

SAMPLE

Serum

a) Collection: obtain serum in the usual way. Do not use plasma. Prostatic acid phosphatase is extremely unstable "in vitro" and at room temperature, due to serum pH, enzymatic activity could decrease up to 50% in a few hours. Therefore, serum should be separated from clot within one hour from collection and stored at 2-10°C.

b) Additives: to avoid inactivation during storage acidify the sample adding 50 ul Reagent D (0.7 mol/l acetic) per 1 ml of serum.

c) Known interfering substances:

- Avoid the use of sera with evident jaundice because the enzymatic activity has low recovery.
- No interferences have been observed with lipemia up to 600 mg/dl triglycerides, or hemoglobin up to 100 mg/dl.
- Avoid the use of plasma, anticoagulants interfere in the reaction.

Refer to Young, D.S. under references for drugs' effect on the present method.

d) Stability and storage instructions: if the preservative detailed in b) is used, sample may be stored at 2-10°C for several days without significant activity loss. Otherwise, proceed as described in a).

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Micropipettes and pipettes for measuring stated volumes.
- Square spectrophotometric cuvettes.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 409 nm
- Reaction temperature: 37°C
- Reaction time: 3 minutes

PROCEDURE

In a cuvette at 37°C, place:

	ACP	ACP-NP
Reagent B ACP	1 ml	-
Reagent B ACP-NP	-	1 ml

Preincubate for a few minutes. Then add:

Sample	100 ul	100 ul
--------	--------	--------

ACP: mix and simultaneously start stopwatch. At 5 minutes measure O.D. Then read absorbance every minute during 3 minutes. Determine average absorbance/minute ($\Delta A/\text{min}$) subtracting each reading from the previous one and averaging values. Use this mean for calculations ($\Delta A_1/\text{min}$).

ACP-NP: proceed as indicated for ACP. $\Delta A_2/\text{min}$ is obtained.

CALCULATIONS

Total Acid Phosphatase (ACP) = $\Delta A_1/\text{min} \times 743 = (U/l_1)$

Non-Prostatic Acid Phosphatase (ACP-NP) =

$\Delta A_2/\text{min} \times 743 = (U/l_2)$

Prostatic Acid Phosphatase (PAP) (U/l) = $U/l_1 - U/l_2$

REFERENCE VALUES

Total Acid Phosphatase (ACP): < 9 U/l (at 37°C)

Prostatic Acid Phosphatase (PAP): < 3.5 U/l (at 37°C)

ACP activity is rare in healthy men and almost absent in women. Thus, expected values are within the limit of instrument error or very close to the system's detection limit. Nonetheless, it is recommended that each laboratory establish its own reference values.

PROCEDURE LIMITATIONS

- See Known interfering substances in SAMPLE.
- The enzymes are sensitive to the action of certain contaminants and enzymatic poisons (heavy metals, cyanides, surfactants, etc.). Thus, it is recommended to maintain strict precautions when cleaning the material used in sample collection and determination.
- Iatrogenic alterations: some medicines may affect the PAP plasma values, as well as massage, catheterization, and other prostate manipulations. Therefore, it is convenient to ask the patient or the physician about the medication,

treatment or diagnosis procedure to which the patient is subjected.

- Clinical diagnosis should not be performed based on the result of a single assay, but it should incorporate clinical and laboratory data.

QUALITY CONTROL METHOD

Each time the test is performed, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known total and prostatic acid phosphatases activity.

PERFORMANCE

a) Reproducibility: CLSI protocol EP15A was applied. Two activity levels were analyzed, in replicates of four, during 5 days, using autoanalyzer. With the obtained data, total and intra-assay precision were calculated.

Total Acid Phosphatase

Intra-assay precision

Level	S.D.	C.V.
23.7 U/l	± 0.400 U/l	1.7 %
45.2 U/l	± 0.575 U/l	1.3 %

Total precision

Level	S.D.	C.V.
23.7 U/l	± 0.480 U/l	2.0 %
45.2 U/l	± 0.872 U/l	1.9 %

Non-Prostatic Acid Phosphatase

Intra-assay precision

Level	S.D.	C.V.
11.6 U/l	± 0.307 U/l	2.7 %
24.5 U/l	± 0.718 U/l	2.9 %

Total precision

Level	S.D.	C.V.
11.6 U/l	± 0.351 U/l	3.0 %
24.5 U/l	± 0.776 U/l	3.2 %

b) Linearity: reaction is linear up to 74 U/l ($\Delta A/\text{min} = 0.100$ O.D.). For higher values repeat determination with diluted sample (0.1 ml sample + 0.2 ml saline solution), multiplying the obtained result by 3.

c) Detection limit: the minimum detectable change of total acid phosphatase activity will be 0.51 U/l, and 0.48 U/l for non-prostatic acid phosphatase.

WIENER LAB. PROVIDES

- 1 x 40 ml Reagent A
- 20 x 2 ml Reagent B
- 1 x 40 ml Reagent C
- 1 x 5 ml Reagent D (Cat. 1351402)

- 1 x 40 ml Reagent A
- 20 x 2 ml Reagent B
- 1 x 40 ml Reagent C
- 1 x 5 ml Reagent D (Cat. 1009669)

- 1 x 40 ml Reagent A
 - 20 x 2 ml Reagent B
 - 1 x 40 ml Reagent C
 - 1 x 5 ml Reagent D
- (Cat. 1009922)*

REFERENCES

- Hillmann, G.Z. - Klin. Chem. Klin. Biochem. 9:273 (1971).
- Skelley, D.S. - Ligand Rev. 2/1:43 (1980).
- Shaw, L.M.; Brummund, W.; Dorio, R.J. - Am. J. Clin. Pathol. 68/1:57 (1977).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., (2001)
- CLSI: Clinical and Laboratory Standards Institute (ex-NCCLS) - Protocols EP 15A, 2001 / EP 17A, (2004).
- Junge, W.; Thormeyer, I.; Schlottmann A. et al - 3rd Alpe-Adria Congress on Clinical Chemistry and Laboratory Medicine, September 7-9 (1994).

SYMBOLS

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device

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
(Non-Prostatic Acid Phosphatase, ACP-NP)

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Wiener Laboratorios S.A.I.C.