# Lipasa

## AA

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Kinetic method for lipase determination in serum and plasma

## SUMMARY

The lipase is mainly produced by the exocrine pancreas and is also being secreted in smaller quantities by the salivary glands as well as by gastric, intestinal and pulmonary mucous membranes. Lipases hydrolyse glycerol esters of fatty acids. Lipase determination is useful in the diagnosis and treatment of pancreatic pathologies such as acute pancreatitis and pancreatic duct's obstruction. Clinical diagnosis should be performed considering the patient's anamnesis, clinical records and additional laboratory results.

## PRINCIPLE

Lipase hydrolyses the specific substrate 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin) ester to release glutaric acid-(6'-methylresorufin) ester, which is an unstable compound spontaneously decomposed releasing a colored compound (methylresorufin) measured at 570 nm. The color intensity is directly proportional to the enzymatic activity.

## PROVIDED REAGENTS

**A. Reagent A:** 50 mmol/l Goods buffer, pH 8.0 with colipase and bile salts.

**B. Reagent B:** 1,2-O-dilauryl-rac-glycerol-3-glutaric acid-(6'-methylresorufin) ester (lipase substrate) in 10 mmol/l tartrate buffer.

## NON-PROVIDED REAGENTS

- Wiener lab's **Calibrador A plus**.
- Saline solution (9 g/l NaCl).

## INSTRUCTIONS FOR USE

**Provided Reagents:** ready to use.

## WARNING

Reagents are for "in vitro" diagnostic use. Use the reagents with the normal precautions for handling all laboratory materials. All the reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents are stable at 2-10°C until the expiration date stated on the box.

Protect from direct light.

## INSTABILITY OR DETERIORATION OF REAGENTS

Reagent A is limpid fluid. Discard in case of turbidity.

Reagent B is a slightly opalescent orange microsuspension. Discard in the presence of a clearly red color.

## SAMPLE

Serum or heparinized plasma

**a) Collection:** obtain serum in the usual way. Separate immediately serum from clot.

**b) Additives:** in case plasma is used as sample, use heparin for its collection.

**c) Known interfering substances:** no interference is observed by bilirubin up to 40 mg/dl, hemoglobin up to 500 mg/dl and triglycerides up to 1200 mg/dl.

Refer to Young bibliography for effects of drugs on this method.

**d) Stability and storage instructions:** lipase in serum or plasma is stable for up to one week at 2-10°C and for up to one year at -20°C.

## REQUIRED MATERIAL (non-provided)

- Volumetric material for measuring stated volumes.
- Autoanalyzer.

## PROCEDURE

(Autoanalyzers)

The general testing procedure for **Lipasa AA líquida** in an autoanalyzer is detailed below. For programming instructions check the user's manual of the autoanalyzer in use.

<b>Sample or Calibrator</b>	2 ul
-----------------------------	------

<b>Reagent A</b>	100 ul
------------------	--------

Incubate for 300 seconds at 37°C.

<b>Reagent B</b>	25 ul
------------------	-------

Incubate for 90 seconds at 37°C. Measure initial absorbance at 575 nm ( $A_1$ ). At exactly 60 seconds measured (using stopwatch) the absorbance again ( $A_2$ ). To obtain lipase result in U/l, multiply the absorbance difference ( $\Delta A = A_2 - A_1$ ) by the factor.

## CALIBRATION

**Calibrador A plus** is processed as a sample and the corresponding factor is calculated based on it. Insert the calibrator's concentration value every time a lot is changed.

## QUALITY CONTROL METHOD

For each determination, process 2 levels of a quality control material (**Standatrol S-E 2 niveles**) with known lipase activities.

## REFERENCE VALUES

Adults: 13 - 60 U/l (37°C)

Each laboratory should establish its own reference values.

## SI SYSTEM UNITS CONVERSION

Lipase (U/l) x 0.017 = Lipase (ukat/l)

## PROCEDURE LIMITATIONS

See Known Interfering Substances under SAMPLE.

Do not pipette by mouth.

It is recommended to use Wiener lab's **Standatrol S-E 2 niveles** as quality control material. The use of controls from other manufacturers may yield different values for certain ranges because they are dependent from the method or system used.

It is important to avoid carry-over contamination in cuvettes and needles when they have been used for triglyceride, cholesterol and HDL-cholesterol and LDL-cholesterol determinations. Therefore, the use of additional cleaning programs is recommended. It is convenient to perform lipase determination independently from other assays.

## PERFORMANCE

**a) Reproducibility:** the following coefficients of variations were obtained as an estimation of the intra-assay (CV<sub>i</sub>) and total (CV<sub>t</sub>) precision studies, according to the guidelines contained in CLSI EP15-A document:

Level	C.V <sub>i</sub>	C.V <sub>t</sub>
34.6 U/l	± 0.99%	4.0 %
59.7 U/l	± 2.80%	3.85 %
97.0 U/l	± 1.51%	4.03 %

**b) Detection limit:** 2 U/l

**c) Linearity:** reaction is linear up to a concentration of 300 U/l. For higher values, dilute sample 1:10 with saline solution (9 g/l ClNa) and repeat the test multiplying the obtained result by 10.

## PARAMETERS FOR AUTOMATIC ANALYZERS

For programming instructions check the user's manual of the autoanalyzer in use. For calibration, Wiener lab's **Calibrador A plus** should be used.

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Wiener Laboratorios S.A.I.C.