# Prealbumin

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# Prealbumin

Immunoturbidimetric method for quantitative determination of prealbumin in serum

## SUMMARY

Prealbumin, also referred as transthyretin, is a tryptophan-rich protein, synthesized by the liver.

It binds to retinol-binding protein in order to transport vitamin A from the liver to the peripheral tissues. It also functions as thyroxine transport protein.

Prealbumin has a half-life in plasma of about 2 days and it is catabolized in the kidney.

Due to the short half-life, its concentration diminishes rapidly when decreasing hepatic synthesis as a result of liver damage or inadequate intake.

Thus, it is useful as a sensitive indicator of malnutrition and of early response to enteral or parenteral nutritional treatment. Prealbumin levels reflect recent dietary intake rather than overall nutritional status.

Because of its synthesis in the liver prealbumin is a reliable index of liver function in hepatobiliary diseases.

This protein is increased in uremia and dehydration, and decreased by an inflammatory response, fasting, hyperthyroidism, severe liver injury, over-hydration and pregnancy.

## PRINCIPLE

Prealbumin reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to prealbumin concentration in the sample and may be spectrophotometrically measured.

## PROVIDED REAGENTS

**A. Reagent A:** phosphate buffer, pH 7.4.

**B. Reagent B:** polyclonal antibodies anti-human prealbumin (goat) in phosphate buffer, pH 7.4.

## NON-PROVIDED REAGENTS

- Saline solution
- Wiener lab.'s **Calibrador Proteínas nivel alto Turbitest AA**

## INSTRUCTIONS FOR USE

**Provided Reagents:** ready to use.

## WARNINGS

The reagents are for "in vitro" diagnostic use.

All patient samples should be handled as though capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable at 2-10°C until the expiration date stated on the box. Do not freeze.

## SAMPLE

Serum

**a) Collection:** obtain in the usual way.

**b) Additives:** not required.

**c) Known interfering substances:** do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation.

No interferences have been observed with hemoglobin up to 1000 mg/dl, triglycerides up to 2500 mg/dl, bilirubin up to 20 mg/dl, heparin up to 50 mg/dl and sodium citrate up to 1000 mg/dl.

See Young, D.S. in References for effect of drugs on the present method.

**d) Stability and storage instructions:** sample should be preferably fresh. In case it cannot be processed immediately, sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

## REQUIRED MATERIAL (non-provided)

- Spectrophotometer
- Square spectrophotometric cuvettes
- Micropipettes and pipettes for measuring the stated volumes
- Kahn or hemolysis tubes
- Stopwatch

## ASSAY CONDITIONS

- Wavelength: 340 nm
- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.
- Reaction time: 15 minutes
- Sample volume: 10 ul
- Final reaction volume: 1810 ul

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

## PROCEDURE

### CALIBRATION CURVE

In Kahn tubes dilute the Calibrador Proteínas nivel alto with saline solution 1:1, 1:2, 1:4, 1:8 and 1:16, using saline solution as the zero point.

<b>Diluted Calibrador Proteínas</b>	10 ul
<b>Reagent A</b>	1500 ul

Homogenize and measure the absorbance of each dilution at 340 nm (OD<sub>1</sub>), setting the instrument to zero with distilled water. Then, add:

<b>Reagent B</b>	300 ul
------------------	--------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD<sub>2</sub>), setting the instrument to zero with distilled water.

Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each Calibrador Proteínas dilution, including the zero point. Draw on graph paper the absorbance differences ( $\Delta A$ ) based on the Calibrador Proteínas concentration in mg/dl (g/l).

#### SAMPLES PROCEDURE

<b>Sample</b>	10 ul
---------------	-------

<b>Reagent A</b>	1500 ul
------------------	---------

Homogenize and measure the absorbance at 340 nm (OD<sub>1</sub>), setting the instrument to zero with distilled water. Then add:

<b>Reagent B</b>	300 ul
------------------	--------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD<sub>2</sub>), setting the instrument to zero with distilled water.

prealbumin levels were assayed and the following results were obtained:

#### Intra-assay precision

Level	S.D.	C.V.
7.7 mg/dl	± 0.3 mg/dl	4.2%
31.9 mg/dl	± 0.7 mg/dl	2.1%
71.5 mg/dl	± 1.7 mg/dl	2.4%

**b) Detection limit:** 4 mg/dl.

**c) Measuring range:** 4 - 80 mg/dl.

**d) Prozone effect:** not noted until 150 mg/dl prealbumin.

#### WIENER LAB. PROVIDES

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009360)

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009658)

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009967)\*

#### REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.

- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5<sup>th</sup> ed., 2000.

#### CALCULATIONS

Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each sample tested. Interpolate this  $\Delta A$  in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the Calibrador Proteínas nivel alto must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

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The Control should be processed in the same manner as the samples.

#### REFERENCE VALUES

20 - 40 mg/dl (0.2- 0.4 g/l)

Each laboratory should set its own reference values.

#### PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control. Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

#### PERFORMANCE

**a) Reproducibility:** replicates of samples containing different

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Wiener Laboratorios S.A.I.C.