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- Expert Panel of National Cholesterol Education Program - JAMA 285/19:2486 (2001).

- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

TG Color

GPO/PAP AA

Enzymatic method for the determination of triglycerides
in serum or plasma




SUMMARY

Triglycerides are lipids absorbed from the diet and also those endogenously produced from carbohydrates.

Their measurement is important for the diagnosis and control of hyperlipemia. These diseases may have a genetic origin or be secondary to others such as nephrosis, diabetes mellitus and endocrinous dysfunctions. The increase of triglycerides has been identified as a risk factor for atherosclerotic diseases.

PRINCIPLE

The reaction system is as follows:

PROVIDED REAGENTS

A. Reagent A: solution containing Good buffer (pH 6.8), chlorophenol, lipoprotein lipase (LPL), glycerol kinase (GK), glycerol phosphate oxidase (GPO), peroxidase (POD), adenosine triphosphate (ATP) and 4-aminophenazone (4-AP).
S. Standard*: 2.26 mmol/l glycerol solution (it is equivalent to 2 g/l triolein).

Final concentrations

Good.....	50 mmol/l, pH 6.8
Chlorophenol.....	2 mmol/l
Lipoprotein lipase.....	≥ 800 U/l
GK.....	≥ 500 U/l
GPO.....	≥ 1500 U/l
POD.....	≥ 900 U/l
ATP.....	2 mol/l
4-AP.....	0.4 mmol/l

NON-PROVIDED REAGENTS

Wiener lab.'s **Calibrador A plus**.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.
Use the reagents according to the working procedures for

clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator (2-10°C) and protected from light until the expiration date stated on the box. Do not expose to high temperatures for long periods of time.

INSTABILITY OR DETERIORATION OF REAGENTS

The Reagent A may develop a slight pink coloration that does not affect its performance.

Blank readings over 0.250 O.D. or abnormally low Standard readings indicate Reagent deterioration. Discard in such case.

SAMPLE

Serum or plasma

a) Collection: obtain serum or plasma after 12-14 hours fasting. Separate red blood cells within 2 hours after extraction.

b) Additives: it is recommended to use Wiener lab's **Anticoagulante W** or heparin as anticoagulant to obtain it.

c) Known interfering substances: no interference is observed from bilirubin up to 15 mg/dl, nor intense hemolysis. See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: triglycerides in serum are stable for 3 days in refrigerator (2-10°C). Do not freeze.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Micropipettes and pipettes for measuring the stated volumes
- Spectrophotometric cuvettes.
- Water bath at 37°C
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 505 nm
- Reaction temperature: 37°C
- Reaction time: 5 minutes
- Sample volume: 10 ul
- Reagent volume: 1 ml
- Final reaction volume: 1.01 ml

PROCEDURE

Homogenize the sample before using, particularly when

dealing with turbid sera.

In three spectrophotometer cuvettes, labeled B (Blank), S (Standard) and U (Unknown), place:

	B	S	U
Sample	-	-	10 ul
Standard	-	10 ul	
Reagent A	1 ml	1 ml	1 ml

Mix. Incubate 5 minutes at 37°C, or 20 minutes at room temperature (18-25°C). Let cool and read in spectrophotometer at 505 nm, setting the instrument to zero with distilled water.

Microtechnique

Follow the previously outlined procedure but using 5 ul Sample and 500 ul Reagent A.

STABILITY OF FINAL REACTION

Final reaction color is stable for 60 minutes, thus absorbance should be read within that period.

CALCULATIONS

Correct readings with the Reagents Blank and use these readings for calculations.

$$\text{TG g/l} = \text{U} \times \text{factor} \quad \text{f} = \frac{2 \text{ g/l}}{\text{S}}$$

QUALITY CONTROL METHOD

Each time the test is performed, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known triglycerides concentration.

UNITS CONVERSION

Triglycerides (g/l) = 0.01 x Triglycerides (mg/dl)

Triglycerides (mg/dl) x 0.0113 = Triglycerides (mmol/l)

REFERENCE VALUES

The National Cholesterol Education Program (NCEP) expert panel provides the following Triglycerides values:

Recommended: < 1.50 g/l

Slightly high to high: 1.50-1.99 g/l

High: 2.00-4.99 g/l

Very high: ≥ 5.00 g/l

However, each laboratory should establish its own references values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Reducing agents decrease the color response, while oxidants color the Reagent increasing the Blanks.

Contaminations with glycerol produce falsely increased results.

PERFORMANCE

a) Reproducibility: processing 20 replicates of the same

samples simultaneously, the following values were obtained:

Level	S.D.	C.V.
0.76 g/l	± 0.039 g/l	0.50 %
3.73 g/l	± 0.060 g/l	0.16 %

b) Recovery: adding known amounts of triolein to different sera, for the whole linearity range of the method, a recovery between 96 and 101% was obtained.

c) Linearity: the reaction is linear up to 10 g/l triglycerides. For higher values, dilute 1:2 with saline solution and repeat the test. Multiply the obtained result by the dilution performed.

d) Detection limit: it depends on the photometer used. In spectrophotometers, for a 0.001 O.D. absorbance variation, the minimum concentration variation detected for the given assay conditions, will be of approximately 0.009 g/l.

PARAMETERS FOR AUTOANALYZERS

For programming instructions check the user's manual of the autoanalyzer in use. For calibration use Wiener lab's **Calibrador A plus**, following the autoanalyzer requirements.

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Wiener Laboratorios S.A.I.C.