# $\alpha$ 2-Macroglobulin

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60 ml: - 1 x 50 ml Reagente A  
 - 1 x 10 ml Reagente B  
 (Cód. 1009945)\*

#### REFERÊNCIAS

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACCC Press, 5<sup>th</sup> ed., 2000.
- Tietz Textbook of Clinical Chemistry - Burtis, C.; Ashwood, E. (5<sup>th</sup> Edition) WB Saunders, 2001.

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# $\alpha$ 2-Macroglobulin

Immunoturbidimetric method for quantitative determination of  $\alpha$ 2-macroglobulin in serum

## SUMMARY

$\alpha$ 2-macroglobulin (AMG) is a glycoprotein, produced by liver which belongs to the group of protease inhibitor proteins. AMG inhibits many different classes of proteases including serine, cysteine and metallo-proteases.

AMG has an important role in the inactivation of plasmin after  $\alpha$ 2-antiplasmin consumption.

Markedly increased AMG levels may be found in individuals with nephrotic syndrome.

Decreased AMG levels are found in individuals with severe acute pancreatitis, liver tumors, disseminated intravascular coagulation, malnutrition and gastroenteropathies.

## PRINCIPLE

The  $\alpha$ 2-macroglobulin reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to the  $\alpha$ 2-macroglobulin concentration in sample and may be spectrophotometrically measured.

## PROVIDED REAGENTS

**A. Reagent A:** phosphate buffer, pH 7.4.

**B. Reagent B:** polyclonal antibodies anti-human AMG (goat) in phosphate buffer, pH 7.4.

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**Provided Reagents:** ready to use.

## WARNINGS

Reagents are for "in vitro" diagnostic use.

All patient samples should be handled as though capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable at 2-10°C until the expiration date stated on the box. Do not freeze.

## SAMPLE

Serum or heparinized plasma

**a) Collection:** obtain in the usual way.

**b) Additives:** when using plasma, heparin is recommended as anticoagulant.

**c) Known interfering substances:** do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation. No interferences are observed with hemoglobin up to 1100 mg/dl, triglycerides up to 2500 mg/dl, bilirubin up to 13 mg/dl and rheumatoid factor up to 300 IU/ml.

See Young, D.S. in References for effect of drugs on the present method.

**d) Stability and storage instructions:** sample should be preferably fresh. In case it cannot be processed immediately, sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

## REQUIRED MATERIAL (non-provided)

- Spectrophotometer.

- Square spectrophotometric cuvettes.

- Micropipettes and pipettes for measuring the stated volumes.

- Kahn or hemolysis tubes.

- Stopwatch.

## ASSAY CONDITIONS

- Wavelength: 340 nm

- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.

- Reaction time: 15 minutes

- Sample volume: 70  $\mu$ l

- Final reaction volume: 1870  $\mu$ l

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

## PROCEDURE

### CALIBRATION CURVE

In Kahn tubes dilute the **Calibrador Proteínas nivel alto** with saline solution 1:10, 1:20, 1:40, 1:80 and 1:160, using saline solution as the zero point.

<b>Diluted Calibrador Proteínas</b>	70 $\mu$ l
<b>Reagent A</b>	1500 $\mu$ l
Homogenize and measure the absorbance of each dilution at 340 nm (OD <sub>1</sub> ), setting the instrument to zero with distilled water. Then, add:	
<b>Reagent B</b>	300 $\mu$ l

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm ( $OD_2$ ), setting the instrument to zero with distilled water. Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each Calibrador Proteínas dilution, including the zero point.

Draw on graph paper the  $\Delta A$  absorbance differences based on the Calibrador Proteínas concentration in mg/dl (g/l).

#### SAMPLES PROCEDURE

Dilute the Samples 1:10 with saline solution.

<b>Diluted Sample</b>	70 $\mu$ l
-----------------------	------------

<b>Reagent A</b>	1500 $\mu$ l
------------------	--------------

Homogenize and measure the absorbance at 340 nm ( $OD_1$ ), setting the instrument to zero with distilled water. Then add:

<b>Reagent B</b>	300 $\mu$ l
------------------	-------------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm ( $OD_2$ ), setting the instrument to zero with distilled water.

#### Intra-assay precision

Level	S.D.	C.V.
44.3 mg/dl	$\pm 2.1$ mg/dl	4.8 %
173.5 mg/dl	$\pm 3.3$ mg/dl	1.9 %
401.8 mg/dl	$\pm 21.3$ mg/dl	5.3 %

**b) Detection limit:** 10 mg/dl.

**c) Measuring range:** 10 - 460 mg/dl.

**d) Prozone effect:** not noted until 4660 mg/dl AMG.

#### WIENER LAB. PROVIDES

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009352)

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009945)\*

#### REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5<sup>th</sup> ed., 2000.
- Tietz Textbook of Clinical Chemistry - Burtis, C.; Ashwood, E. (5<sup>o</sup> Edition) WB Saunders, 2001.

#### CALCULATIONS

Calculate the absorbance difference ( $\Delta A = OD_2 - OD_1$ ) for each sample tested. Interpolate this  $\Delta A$  in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the **Calibrador Proteínas nivel alto** must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

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Control should be processed in the same manner as samples.

#### REFERENCE VALUES

130 - 300 mg/dl (1.3 - 3.0 g/l)

Each laboratory should set its own reference values.

#### PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control.

Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

#### PERFORMANCE

**a) Reproducibility:** replicates of samples containing different AMG levels were assayed and the following results were obtained:

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Wiener Laboratorios S.A.I.C.