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60 ml: - 1 x 50 ml Reagente A
- 1 x 10 ml Reagente B
(Cód. 1009948)*

REFERÊNCIAS

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.

Apo B

Immunoturbidimetric method for quantitative determination of apolipoprotein B in serum

SUMMARY

Apolipoprotein B (Apo B) is one of the larger human proteins. Apo B represents approximately 95% of the total protein content of low density lipoprotein (LDL). It is also the protein component of very low density lipoprotein (VLDL) and chylomicrons.

Apo B is essential in lipoproteins formation and releases them to plasma. Apo B interacts with LDL receptors present in peripheral cells resulting in removal and degradation of LDL. Apo B elevated levels are found in patients with familial hypercholesterolemia, familial hiperapobetalipoproteinemia, nephrotic syndrome, biliary obstruction, type II hyperlipidemia. Increased Apo B levels are associated with increased risk of atherosclerosis.

Decreased Apo B levels are found in patients with familial hypobetalipoproteinemia, abetalipoproteinemia, chronic anemia, liver disease and neuromuscular degeneration.

PRINCIPLE

The apolipoprotein B reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to the apolipoprotein B concentration in the sample and may be spectrophotometrically measured.

PROVIDED REAGENTS

A. Reagent A: phosphate buffer, pH 7.4.

B. Reagent B: polyclonal antibodies anti-human apolipoprotein B (goat) in phosphate buffer, pH 7.4.

NON-PROVIDED REAGENTS

- Saline solution.

- Wiener lab.'s **Apo Calibrator Turbitest AA.**

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

The reagents are for "in vitro" diagnostic use.

All patient samples should be handled as though capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box. Do not freeze.

SAMPLE

Serum

a) Collection: obtain in the usual way.

b) Additives: not required.

c) Known interfering substances: do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation. No interferences are observed with hemoglobin up to 1000 mg/dl, triglycerides up to 2500 mg/dl, bilirubin up to 20 mg/dl, heparin up to 50 mg/dl and sodium citrate up to 1000 mg/dl.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: sample should be preferably fresh. In case it cannot be processed immediately, the sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.

- Square spectrophotometric cuvettes.

- Micropipettes and pipettes for measuring the stated volumes.

- Kahn or hemolysis tubes.

- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 340 nm

- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.

- Reaction time: 15 minutes

- Sample volume: 10 ul

- Final reaction volume: 1810 ul

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

PROCEDURE

CALIBRATION CURVE

In Kahn tubes dilute the Apo Calibrator Turbitest AA with saline solution 1:1, 1:2, 1:4, 1:8 and 1:16, using saline solution as the zero point.

Diluted Apo Calibrator	10 ul
------------------------	-------

Reagent A	1500 ul
-----------	---------

Homogenize and measure the absorbance of each dilution at 340 nm (OD₃₄₀), setting the instrument to zero with distilled water. Then, add:

Reagent B	300 ul
Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD ₂), setting the instrument to zero with distilled water. Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each Apo Calibrator dilution, including the zero point. Draw on graph paper the ΔA absorbance differences based on the Apo Calibrator concentration in mg/dl (g/l).	
SAMPLES PROCEDURE	
Sample	10 ul
Reagent A	1500 ul
Homogenize and measure the absorbance at 340 nm (OD ₁), setting the instrument to zero with distilled water. Then add:	
Reagent B	300 ul
Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD ₂), setting the instrument to zero with distilled water.	

Intra-assay precision

Level	S.D.	C.V.
133.7 mg/dl	± 3.7 mg/dl	2.8%
169.1 mg/dl	± 4.2 mg/dl	2.5%

b) Detection limit: 30 mg/dl

c) Measuring range: 30 - 400 mg/dl

d) Prozone effect: not noted until 4000 mg/dl Apo B.

WIENER LAB. PROVIDES

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009362)

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009640)

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009948)*

REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.

CALCULATIONS

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each sample tested. Interpolate this ΔA in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the Apo Calibrator Turbitest AA must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

QUALITY CONTROL METHOD

It is recommended the use of **Lipid Control**, separately provided by Wiener lab. The control should be processed in the same manner as samples.

REFERENCE VALUES

Men: 79 - 127 mg/dl (0.79 - 1.27 g/l)

Women: 70 - 122 mg/dl (0.70 - 1.22 g/l)

Each laboratory should set its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE. It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control. Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

PERFORMANCE

a) Reproducibility: replicates of samples containing different apolipoprotein B levels were assayed and the following results were obtained:

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Wiener Laboratorios S.A.I.C.