Cystatin C

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Cystatin C

Immunoturbidimetric method for the determination of cystatin C

SUMMARY

Cystatin C is a low molecular weight (13 kDa) cytoplasmic protein, functioning as an inhibitor of various cystatin proteases. Cystatin C has a stable production rate and is removed from the bloodstream by glomerular filtration.

In healthy individuals cystatin C is completely reabsorbed and degraded in the tubules but in subjects with renal disorders its level in blood may be raised as high as 2 to 5 times normal values. Unlike creatinine, cystatin C is unaffected by inflammatory processes, sex, age, diet, and nutritional status. Numerous studies have shown that serum cystatin C is a better marker of Glomerular Filtration Rate (GFR) than serum creatinine.

Cystatin C measurements are useful in the diagnosis and treatment of renal diseases.

PRINCIPLE

Cystatin C reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to the cystatin C concentration in the sample and may be spectrophotometrically measured.

PROVIDED REAGENTS

A. Reactivo A: Tris buffer 100 mM, pH 8.5 ± 0.3.

B. Reactivo B: synthesized polystyrene latex particles coated with polyclonal anti-cystatin C antibodies (rabbit) in glycine buffer 100 mM, pH 7.3 ± 0.3.

NON PROVIDED REAGENTS

- Wiener lab's **Cystatin C Calibrator Turbitest AA**.
- NaCl 9% saline solution.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

Reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box.

Once opened, reagents are stable for up to 30 days store at 2-10°C. Do not freeze.

SAMPLE

Serum or plasma

a) Collection: obtain serum in the usual way. Separate immediately serum from clot.

b) Additives: in case plasma is used as sample, employ heparin or EDTA for its collection.

c) Known interfering substances: hemoglobin up to 500 mg/dl, bilirubin up to 30 mg/dl, triglycerides up to 5 g/dl and ascorbic acid up to 50 mg/dl do not interfere with cystatin C determination.

See Young, S.D. in references for effect of drugs on the present method.

d) Stability and storage instructions: after sampling, the test should be performed without delay. If the test cannot be performed immediately, sample should be placed in a tightly sealed container and store for up to 1 week at 2-10°C, for up to 2 days at 20-25°C or for up to three months at -20°C. Do not repeatedly freezing and thawing.

MATERIAL REQUIRED (non-provided)

- Volumetric material for measuring stated volumes.
- Automated analyzer.

PROCEDURE

(Automated analyzer)

The general testing procedure for cystatin C in an automated analyzer is detailed below.

For programming instructions check the user's manual of the automated analyzer in use.

Sample or Calibrator	3 ul
Reagent A	230 ul
Incubate during 300 segundos a 37°C.	
Reagent B	50 ul

Incubate during 4 minutes, 30 seconds and then measure absorbance change at 546 nm.

Calculate cystatin C concentration with the measure absorbance change by interpolation from a calibration curve prepared with calibrators of known concentrations.

CALIBRATION

Cystatin C Calibrator Turbitest AA provided separately by Wiener lab. should be used for calibration and 0.9% NaCl shall be used for zero calibrator (blank solution).

It is recommended that each laboratory determine calibration

frequency, as this would depend on the analyzer in use as well as the types and number of other assays being run. A new calibration curve should be drawn at least once a month or when a new lot of reagent is used.

QUALITY CONTROL METHOD

For each assay, process 2 levels of a quality control material (**Cystatin C Control Turbitest AA**) with known cystatin C concentrations

QC intervals and limits should be adapted to each laboratory's individual requirements. Each laboratory should establish corrective measures if values fall outside the limits.

REFERENCE VALUE

Adults: 0.56 - 1.02 mg/l

It is recommended that each laboratory establish its own reference intervals.

PROCEDURE LIMITATIONS

See Known Interfering Substances under SAMPLE.

Do not pipette by mouth.

Avoid contamination to preserve the integrity of the reagents.

Only use thoroughly clean and dry micropipettes for measurements.

Samples containing cystatin C over the assay range should be diluted with 0.9% NaCl and retested. Correct the results for the dilution factor.

PERFORMANCE

a) Reproducibility: two different levels (high and low) were assayed twenty times in order to obtain intra-assay precision:

Level	S.D.	C.V.
0.84 mg/l	± 0.005 mg/l	0.58 %
4.02 mg/l	± 0.022 mg/l	0.55 %

b) Linearity: the linearity range of the assay is 0.1 - 10 mg/l cystatin C.

c) Analytical sensitivity ≤ 0.1 mg/l cystatin C.

d) Prozone effect: not note until 60 mg/l cystatin C.

PARAMETERS FOR AUTOMATIC ANALYZERS

For programming instructions check the user's manual of the automated analyzer in use. For calibration, Wiener lab's **Cystatin C Calibrator Turbitest AA** should be used.

WIENER LAB. PROVIDES

50 ml: - 1 x 40 ml Reagent A
- 1 x 10 ml Reagent B
(Cat. N° 1009366)

25 ml: - 1 x 20 ml Reagent A
- 1 x 5 ml Reagent B
(Cat. N° 1009646)

25 ml: - 1 x 20 ml Reagent A
- 1 x 5 ml Reagent B
(Cat. N° 1009955)*

REFERENCES

- Filler G, Bökenkamp A, Hofmann W, Le Bricon T, Martínéz-Brú C, Grubb A. Cystatin C as a marker of GFR - history, indications, and future research. Clin. Biochem. 38: 1, 2005.
- Dharnidharka VR, Kwon C, Stevens G. Serum cystatin C is superior to serum creatinine as a marker of kidney function: a meta-analysis. Am J Kidney Dis. 40, 221, 2002.

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Wiener Laboratorios S.A.I.C.