WL 19 E-Z Cleanser

AA

Para el lavado en analizadores hematológicos diferenciales
Para a lavagem em analisadores hematológicos diferenciais

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WL 19 E-Z Cleanser

AA

Cleaning solution for use in differential hematology analyzers

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Nr kat. 1474513

USES

WL 19 E-Z Cleanser AA is an enzymatic cleaning solution for use in differential hematology analyzers: Wiener lab. Counter 19 and Wiener lab. Counter 19 CP.

PROVIDED REAGENTS

WL 19 E-Z Cleanser AA: 3.0-10 g/l proteolytic enzyme, 3.0-5.0 g/l sodium chloride and 1.0-4.0 g/l phosphate buffer solution.

INSTRUCTIONS FOR USE

For it use, see the analyzer operation manual instructions.

WARNINGS

The Reagent is for "in vitro" diagnostic use.
Do not inhale nor ingest. In case of inhaling or ingestion, seek medical advice.
Avoid contact with the eyes and skin.
Use the reagents according to the working procedures for clinical laboratories.
All reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

WL 19 E-Z Cleanser AA: stable at 2-10°C until the expiration date stated on the box.
Over a period of time, the reagent may develop turbidity or color change, not affecting its performance.
Keep container closed when not in use. Discard opened container after 60 days.

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Wiener Laboratorios S.A.I.C.