

WL 19 Con

AA

LOT

Rev.: 004

WL 19 Con AA 3 niveles

PARAMETRO	WIENER LAB			WIENER LAB			WIENER LAB		
	COUNTER 19			COUNTER 19 CP			COUNTER XS 20		
	NIVEL 1	NIVEL 2	NIVEL 3	NIVEL 1	NIVEL 2	NIVEL 3	NIVEL 1	NIVEL 2	NIVEL 3
WBC (x 10 ⁹ /L)	2,2 ± 0,5	8,2 ± 1,0	21,0 ± 2,5	2,2 ± 0,5	8,3 ± 1,0	21,5 ± 2,5	2,2 ± 0,5	8,0 ± 1,0	21,3 ± 2,5
LYMPH # (x 10 ⁹ /L)	1,2 ± 0,3	2,5 ± 0,7	2,7 ± 1,3	1,2 ± 0,3	2,4 ± 0,7	2,8 ± 1,3	1,4 ± 0,3	2,6 ± 0,7	2,9 ± 1,7
MID # (x 10 ⁹ /L)	0,2 ± 0,2	0,7 ± 0,7	1,2 ± 1,1	0,2 ± 0,2	0,7 ± 0,7	1,1 ± 1,0	0,1 ± 0,1	0,4 ± 0,4	1,2 ± 0,9
GRAN # (x 10 ⁹ /L)	0,7 ± 0,3	5,1 ± 0,7	17,1 ± 1,7	0,8 ± 0,3	5,2 ± 0,7	17,6 ± 1,8	0,7 ± 0,3	5,0 ± 0,7	17,3 ± 2,0
LYMPH % (%)	55,5 ± 12,0	30,0 ± 8,0	13,0 ± 6,0	56,0 ± 12,0	29,0 ± 8,0	13,0 ± 6,0	63,0 ± 12,0	32,2 ± 8,0	13,6 ± 8,0
MID % (%)	11,0 ± 9,0	8,0 ± 7,0	5,5 ± 5,0	9,5 ± 9,0	8,0 ± 7,0	5,0 ± 4,5	5,1 ± 4,0	5,0 ± 4,0	5,4 ± 4,0
GRAN % (%)	33,5 ± 10,0	62,0 ± 8,0	81,5 ± 8,0	34,5 ± 10,0	63,0 ± 8,0	82,0 ± 8,0	31,9 ± 10,0	62,8 ± 8,0	81,0 ± 9,0
HGB (g/dL)	6,2 ± 0,6	13,7 ± 0,9	18,7 ± 1,3	6,1 ± 0,6	13,7 ± 0,9	18,9 ± 1,3	6,2 ± 0,6	13,8 ± 0,9	18,8 ± 1,3
RBC (x 10 ¹² /L)	2,53 ± 0,18	4,82 ± 0,24	5,75 ± 0,30	2,45 ± 0,18	4,76 ± 0,24	5,80 ± 0,30	2,38 ± 0,18	4,55 ± 0,24	5,49 ± 0,30
HCT (%)	19,0 ± 1,5	41,0 ± 2,0	55,2 ± 2,4	18,3 ± 1,5	40,5 ± 2,0	55,7 ± 2,4	17,6 ± 1,5	38,7 ± 2,0	51,9 ± 2,4
MCV (fL)	75,0 ± 5,0	85,0 ± 5,0	96,0 ± 5,0	74,5 ± 5,0	85,0 ± 5,0	96,0 ± 5,0	74,0 ± 5,0	85,0 ± 5,0	94,5 ± 5,0
MCH (pg)	24,5 ± 2,7	28,4 ± 2,5	32,5 ± 2,5	24,9 ± 2,7	28,8 ± 2,5	32,6 ± 2,5	26,1 ± 2,7	30,3 ± 2,5	34,2 ± 2,5
MCHC (g/dL)	32,7 ± 3,0	33,4 ± 3,0	33,9 ± 3,0	33,4 ± 3,0	33,9 ± 3,0	33,9 ± 3,0	35,2 ± 3,0	35,7 ± 3,0	36,2 ± 3,0
RDW-CV (%)	16,0 ± 3,0	15,5 ± 3,0	14,5 ± 3,0	14,5 ± 3,0	14,0 ± 3,0	13,0 ± 3,0	16,5 ± 3,0	15,5 ± 3,0	14,5 ± 3,0
RDW-SD (fL)	36,0 ± 6,0	39,0 ± 6,0	41,0 ± 8,0	36,5 ± 6,0	39,0 ± 6,0	41,0 ± 8,0	43,5 ± 6,0	45,0 ± 6,0	47,0 ± 8,0
PLT (x 10 ⁹ /L)	76 ± 20	268 ± 40	505 ± 60	70 ± 20	258 ± 40	494 ± 60	60 ± 20	266 ± 40	534 ± 60
MPV (fL)	8,6 ± 3,0	8,1 ± 3,0	8,0 ± 3,0	8,0 ± 3,0	8,0 ± 3,0	7,7 ± 3,0	8,5 ± 3,0	9,1 ± 3,0	9,1 ± 3,0

WL 19 Con AA nivel 2

PARAMETRO	NIVEL 2		
	WIENER LAB COUNTER 19	WIENER LAB COUNTER 19 CP	WIENER LAB COUNTER XS 20
WBC (x 10 ⁹ /L)	8,2 ± 1,0	8,3 ± 1,0	8,0 ± 1,0
LYMPH # (x 10 ⁹ /L)	2,5 ± 0,7	2,4 ± 0,7	2,6 ± 0,7
MID # (x 10 ⁹ /L)	0,7 ± 0,7	0,7 ± 0,7	0,4 ± 0,4
GRAN # (x 10 ⁹ /L)	5,1 ± 0,7	5,2 ± 0,7	5,0 ± 0,7
LYMPH % (%)	30,0 ± 8,0	29,0 ± 8,0	32,2 ± 8,0
MID % (%)	8,0 ± 7,0	8,0 ± 7,0	5,0 ± 4,0
GRAN % (%)	62,0 ± 8,0	63,0 ± 8,0	62,8 ± 8,0
HGB (g/dL)	13,7 ± 0,9	13,7 ± 0,9	13,8 ± 0,9
RBC (x 10 ¹² /L)	4,82 ± 0,24	4,76 ± 0,24	4,55 ± 0,24
HCT (%)	41,0 ± 2,0	40,5 ± 2,0	38,7 ± 2,0
MCV (fL)	85,0 ± 5,0	85,0 ± 5,0	85,0 ± 5,0
MCH (pg)	28,4 ± 2,5	28,8 ± 2,5	30,3 ± 2,5
MCHC (g/dL)	33,4 ± 3,0	33,9 ± 3,0	35,7 ± 3,0
RDW-CV (%)	15,5 ± 3,0	14,0 ± 3,0	15,5 ± 3,0
RDW-SD (fL)	39,0 ± 6,0	39,0 ± 6,0	45,0 ± 6,0
PLT (x 10 ⁹ /L)	268 ± 40	258 ± 40	266 ± 40
MPV (fL)	8,1 ± 3,0	8,0 ± 3,0	9,1 ± 3,0

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
fecha	data	date	data
parámetro	parametro	parameter	parametry
nivel	nível	level	poziom

